

Limiting factors: photosynthesis

FizziQ activity · High school · 50 min



Scan to open in FizziQ

High school students vary light, temperature or wavelength and follow CO₂ consumption with an SCD40 sensor. FizziQ Connect logs the data on smartphone or tablet.

Learning objectives

- Design an experimental protocol with identified variables (independent, dependent, controlled)
- Formulate a testable hypothesis about the influence of a factor on photosynthesis
- Perform multiple measurements under controlled conditions
- Calculate and compare CO₂ consumption rates under different conditions
- Draw conclusions based on data and relate them to biological mechanisms

Instruments and sensors

CO₂ meter · O₂ meter (optional) · SCD40 sensor (CO₂, temperature, humidity) · O₂ sensor

Materials

- Smartphone or tablet with FizziQ Connect
- SCD40 sensor (CO₂, temperature, humidity)
- O₂ sensor (optional)
- M5 Stack module
- Sealed transparent chamber
- Fresh leaves
- Light sources of varying intensity and color
- Thermometer and water bath
- Note: the protocol remains adaptable to any comparable gas concentration measurement tool.

Protocol

1. Choose the factor you wish to study: light intensity, light color, temperature, or plant type.
2. Identify the variables: the independent variable (the factor you vary), the dependent variable (CO₂ consumption rate), and the controlled variables (everything else held constant).
3. Prepare the basic setup: connect the SCD40 sensor to the M5 Stack module, connect FizziQ Connect, and verify the readings.
4. Perform a first control measurement: place the plants in the sealed chamber in daylight for 10 minutes. Calculate the CO₂ consumption rate (slope of the CO₂ curve).
5. Modify the chosen factor for the second condition. For example: reduce light intensity (move the lamp farther or use a filter).
6. Perform the second measurement under the same conditions (same duration, same plant volume, same chamber), changing only the tested factor.
7. Repeat for at least a third condition (for example: bright, medium, and dim light; or 5°C, 20°C, 35°C).
8. For each condition, calculate the slope of the CO₂(t) curve during the linear phase. This slope represents the CO₂ consumption rate.
9. Plot a summary graph: CO₂ consumption rate versus the studied factor.
10. Compare your results with your initial hypothesis. Conclude on the influence of the factor and propose explanations based on the underlying biology.

Expected results

Results depend on the factor studied. For light intensity: CO₂ consumption rate increases with intensity up to a plateau (pigment saturation). For light color: red and blue light give the highest rates, green the lowest. For temperature: rate increases up to 25-35°C then drops. For plant type: C4 plants (corn) are faster than C3 plants (spinach) at high light and temperature.

Questions

- Why is photosynthesis more efficient under red and blue light than under green light?
- How does temperature influence enzyme activity, and why is there an optimum?
- What is the light compensation point and why does it matter for plant ecology?
- How would increasing CO₂ concentration affect the photosynthesis rate?
- Why is it important to vary only one factor at a time?
- How do your results relate to the global carbon cycle?

Going further

- Study the effect of monochromatic light using red, blue, and green filters to relate results to the chlorophyll absorption spectrum
- Use aquatic plants (Elodea) and count O₂ bubbles as an alternative measure of photosynthesis rate
- Compare sun leaves and shade leaves from the same tree species
- Investigate the effect of CO₂ concentration by adding baking soda to the chamber
- Study the combined effect of two factors (e.g., light and temperature) and create a 2D response surface